

## Dermatology Sample Collection for Skin, Hair and Nails

<b>Order:</b>	<b>Dermatophyte culture Fungal Cultures</b>	<b>Dermatophyte culture Fungal Cultures</b>	
<b>Specimen:</b>	<b>skin</b>  Vigorously swab area	<b>Hair</b>  <b>Nail</b>	
<b>Collection Container:</b>	Aerobic/Anaerobic transport system (e-swab) <b>Lawson # 78</b>  	Sterile container <b>Lawson # 93</b>  	Sterile red top vacutainer <b>Lawson # 86</b>  

This is not all inclusive refer to Laboratory Services catalog or call Pathology Client Services 8-PATH