

Blood and Bone Marrow Bottle Collection

Order:	Blood Culture	AFB Culture Nocardia Culture	Mold Blood Culture	Fungus Culture
Specimen:	Blood	Blood Bone Marrow	Blood Bone Marrow	Bone Marrow
Collection Instructions:	<ul style="list-style-type: none"> Adults: 5-10 mL Pediatric 2.5-10 mL Infants 0.5- 2mL 			
Collection Container:	<p style="text-align: center;">Blood Culture Vial Workday # 127919</p> 		<p style="text-align: center;">Myco/F Lytic Culture Vial Workday # 433003</p> 	

Order **Mold Blood Culture** for blood stream infections due to filamentous fungi such as Fusarium or Histoplasma.

Mold Blood Culture should **NOT** be ordered to diagnose Candida blood stream infections. Use the **Blood Culture** order for Candida infections.