

Dermatology Sample Collection for Skin, Hair and Nails

Order:	Dermatophyte Screen, Fungus Culture		
Specimen:	Skin (Vigorously swab area)	Hair Nail	
Collection Container:	Aerobic/Anaerobic Transport System (Eswab) Workday # 78 	Sterile Container Workday # 93 	Sterile Red Top Vacutainer Workday # 86 

Candida species will **NOT** be detected by this culture. This culture is intended to diagnose common causes of dermatomycosis caused by dermatophytes, Fusarium, Scopulariopsis, and others.

This is not all inclusive refer to VCU Pathology Catalog or call Pathology Client Services 8-PATH