


# Microbiology Testing for HPV (Human Papillomavirus)\* and Co-collection w/ PAP Testing

*Using one sample both cervical cytology and HPV testing can be performed.*

\*FDA approved for primary HPV testing and reports High-Risk HPV genotypes

**Orderable Test Names in Epic = Pap Test, HPV High Risk PCR Testing (Both orders allow for reflex and co-testing)**

**ACCEPTABLE SPECIMEN:** ThinPrep® Liquid-Based system



(ThinPrep® Can be collected with (in order pictured above): Rovers Cervix-Brush Combi OR Cytobrush/Plastic Spatula OR Broom-like Device).

Workday: #218879

#218878

#218877

**Note: The Rovers Cervix-Brush Combi device should not be used during pregnancy.**

**SOURCE AND BODY SITES:** Sources = Cervical or vaginal

**COLLECTION AND TRANSPORT:** Transport ThinPrep vial to the laboratory at room temperature within 24 hours.

For questions regarding specimen collection and handling, contact **Cytology (Pap)** at 628-2591; **Microbiology (HPV)** at 628-4000.