


Microbiology Testing for HPV (Human Papillomavirus)* and Co-collection w/ PAP Testing, Chlamydia trachomatis, Neisseria gonorrhoeae, & Trichomonas vaginalis

*With one specimen all
testing can be performed.*

*FDA approved for primary HPV testing and reports High-Risk HPV genotypes

Epic Orderable Test Names = Pap Test, HPV High Risk PCR Testing (Both orders allow for co-testing)

ACCEPTABLE SPECIMEN: ThinPrep® Liquid-Based Vial



(ThinPrep® Can be collected with (in order pictured above): Rovers Cervix-Brush Combi **OR** Cytobrush/Plastic Spatula **OR** Broom-like Device).

Workday: #218879

#218878

#218877

Note: The Rovers Cervix-Brush Combi device should not be used during pregnancy.

SOURCE AND BODY SITES: Sources = Cervix or Vagina**

COLLECTION AND TRANSPORT: Transport ThinPrep vial to the laboratory at room temperature within 24 hours.

****For Rectal HPV testing, refer to next page.**

HPV High Risk with HPV 16 and 18, Rectal Collection Device utilized by LabCorp

Test Name:	HPV High Risk with HPV 16 and 18, <u>Rectal</u>
Specimen Types:	Rectal
Testing Laboratory:	LABCORP Lab
Collection Device:	ThinPrep® Vial (collect specimen using a Dacron swab)
Special Handling:	Store/Transport specimens at room temperature until delivery to the laboratory.

Specimen Types	Collection Devices		Handling/Transport
Rectal	ThinPrep® Vial Workday #218872	Dacron (polyester tip) swab Workday #443355	Store/Transport at room temperature.
	