

Mycoplasma genitalium PCR Testing Collection Device

Workday # 441287

Test Name:	Mycoplasma genitalium PCR Testing
Specimen Types:	Endocervix, Vagina
Collection Container:	cobas PCR Media Uni Swab
Special Handling:	Store/Transport specimens at room temperature until delivery to the laboratory.

Specimen Types	Collection Kit Workday # 441287	Handling/Transport
Endocervix Vagina	cobas PCR Media Uni Swab  A photograph of the cobas PCR Media Uni Swab collection kit. It consists of a clear plastic tube with a yellow cap and a blue label, and a white swab with a blue tip. The swab is placed horizontally below the tube.	Store/Transport at room temperature.