






# Microbiology Testing for Sexually Transmitted Disease (STD)

<p><b>Order:</b></p>	<p><i>Chlamydia trachomatis</i>/ <i>Gonococcus (Neisseria gonorrhoeae)</i> by DNA Amplification</p> <p>Trichomonas Direct testing-dedicated tube</p>	<p><i>Chlamydia trachomatis</i>/ <i>Gonococcus (Neisseria gonorrhoeae)</i> by DNA Amplification</p> <p>Trichomonas Direct Testing</p>	<p><i>Gonococcus (Neisseria gonorrhoeae)</i> Culture</p>	<p><i>Chlamydia trachomatis</i> Culture</p>
<p><b>Specimen:</b></p>	<p>Female Endocervix</p>	<p>Female Urine Male Urine First void, 20-50 mL collected into sterile container, using pipette provided in kit add 7 mL of urine into transport tube</p>	<p>Eye Genital Respiratory, Upper</p>	<p>Eye Genital – Penis Genital - Vagina Genital - Vulva Rectal Swab Respiratory, Upper</p>
<p><b>Collection Container:</b></p>	<p>Xpert® Vaginal/ Endocervical specimen collection kit (pink lid) <b>Workday # 129641</b></p> 	<p>Xpert® Urine collection kit (yellow lid) <b>Workday # 39973</b></p> 	<p>Aerobic/Anaerobic transport system (e-swab) <b>Workday # 78</b></p> 	<p>Universal viral transport media (UTM, red lid) <b>Workday # 39659</b></p> 
<p><b>Rectal- refer to next page</b></p>		<p>VCUHealth 12/3/24 AH</p>		

**Chlamydia and Gonococcus, NAA**  
**Collection Device utilized by LabCorp**

**Epic Order:    Miscellaneous Test**

<b>Order Question</b>	
<b>Test code:</b>	188672    (Enter only the reference laboratory test #; Labcorp: 6 numerical digits)
<b>Test Name:</b>	Chlamydia/Gonococcus, <b>Rectal Swab</b> , NAA
<b>Testing Laboratory:</b>	Click LABCORP
<b>Collection Container:</b>	APTIMA® Unisex Swab
<b>Special Handling:</b>	Maintain specimen at room temperature or refrigerate.

Specimen Type	Collection Kit	Handling/Transport
Rectal	<p>APTIMA® Unisex Swab Transport</p> 	Store/Transport at room temperature.