

Microbiology Testing for HSV/VZV, Viral and Chlamydia

Specimens Requiring Universal Transport Media (UTM)

Order:	HSV 1 & 2 Direct testing VZV Direct testing HSV 1 & 2/VZV Direct testing	Virus culture	Chlamydia trachomatis Culture
Specimen:	Source: Lesion (Vesicle) Body sites: Genital Rectal Mouth Any dermal site	Body fluid Eye Herpes Simplex Lower Respiratory (Non-BAL) Stool Tissue Urine Vesicle Varicella Zoster	Eye Genital – Penis Genital - Vagina Genital - Vulva Rectal Swab Respiratory, Upper

Collection Container: Universal viral transport media (UTM, red lid)



Workday # 39659 Contact Material Management

- This chart is not all inclusive, refer to the Pathology Laboratory Services Catalog or call Pathology Client Services & PATH