

# Pap test and HPV Co-collection Labels

Always collect *Pap test* specimen **BEFORE** obtaining specimen for *HPV DNA* testing.

**Order =**

**Pap & HPV**





**Pap w/Reflex HPV**

**Pap (only)**




**HPV (only)**



**Pap = SurePath**



**HPV = STM**



Place extra label in bag with HPV specimen

Bag SurePath and STM specimens separately and include requisitions (if applicable) in both bags

**STM = Specimen Transport Medium™ = hc2 DNA collection device**

Effective October 2014

For questions regarding specimen collection and handling, contact **Cytology (Pap)** at 628-2591 or 628-1288; **Microbiology (HPV)** at 628-4000 or 628-3674.